

Product Datasheet

Euphol (orb1227500)

Description	Euphol
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4°C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	$C_{30}H_{50}O$
Purity	greater than or equal to 98% (HPLC)
MW	426.7
CAS Number	514-47-6
SMILES	<chem>C[C@H](CC/C=C(C)/C)[C@]1([H])CC[C@]([C@@]1(C)CC2)(C)C3=C2[C@]4(C)[C@](CC3)([H])C(C)(C)[C@@H](O)CC4</chem>
Solubility (25°C)	Soluble in DMSO, methanol or ethanol (all 5-10mg/ml). Insoluble in water.
Expiration Date	12 months from date of receipt.