

Product Datasheet

Telaprevir (orb1227503)

Description	Telaprevir
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	C ₃₆ H ₅₃ N ₇ O ₆
Purity	greater than or equal to 95% (HPLC)
MW	679.9
CAS Number	402957-28-2
SMILES	<chem>O=C([C@H](C(C)(C)NC([C@H](C1CCCCC1)NC(C2=NC=CN=C2)=O)=O)N3C(C(N[C@@H](CCC)C(C(NC4CC4)=O)=O)=O)[C@](CCC5)([H])[C@]5([H])C3</chem>
Solubility (25°C)	Soluble in DMSO (10mg/ml) or DMF.
Expiration Date	12 months from date of receipt.

