

Product Datasheet

TDZD-8 (orb1227504)

Description	TDZD-8
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	$C_{10}H_{10}N_2O_2S$
Purity	greater than or equal to 98% (HPLC)
MW	222.3
CAS Number	327036-89-5
SMILES	<chem>O=C(N1CC2=CC=CC=C2)N(C)SC1=O</chem>
Solubility (25°C)	Soluble in DMSO (20mg/ml), DMF (20mg/ml) or ethanol (20mg/ml).
Expiration Date	12 months from date of receipt.

