

## Product Datasheet

### Favipiravir (orb1227532)

<b>Description</b>	Favipiravir
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>5</sub> H <sub>4</sub> FN <sub>3</sub> O <sub>2</sub>
<b>Purity</b>	greater than or equal to 98%
<b>MW</b>	157.1
<b>CAS Number</b>	259793-96-9
<b>SMILES</b>	FC1=CNC(C(C(N)=O)=N1)=O
<b>Solubility (25°C)</b>	Soluble in water (30mg/ml) or DMSO (30mg/ml).
<b>Expiration Date</b>	12 months from date of receipt.

