

Product Datasheet

EIDD-2801 (orb1227535)

Description	EIDD-2801
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4°C. Long Term Storage: -2°C. Handling Advice: Protect from light. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	C ₁₃ H ₁₉ N ₃ O ₇
Purity	greater than or equal to 98% (HPLC)
MW	329.3
CAS Number	2349386-89-4
SMILES	O=C(N(C=C/1)[C@@H]2O[C@H](COC(C(C)C)=O)[C@@H](O)[C@H]2O)N1=N\O
Solubility (25°C)	Soluble in DMSO (20mg/ml), DMF (25mg/ml) or ethanol (25mg/ml).
Expiration Date	12 months from date of receipt.

