

## Product Datasheet

### Hydroxychloroquine . sulfate (orb1227539)

<b>Description</b>	Hydroxychloroquine . sulfate
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{18}H_{26}ClN_3O \cdot H_2SO_4$
<b>Purity</b>	greater than or equal to 98% (HPLC)
<b>MW</b>	335.9 . 98.1
<b>CAS Number</b>	747-36-4
<b>SMILES</b>	<chem>C1C=CC=C(C(NC(CCCN(CC)CCO)C)=CC=N2)C2=C1.O=S(O)(O)=O</chem>
<b>Solubility (25°C)</b>	Soluble in DMSO (30mg/ml) or DMF (30mg/ml).
<b>Expiration Date</b>	12 months from date of receipt.

