

## Product Datasheet

### Camostat . mesylate (orb1227541)

<b>Description</b>	Camostat . mesylate
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Protect from moisture. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{20}H_{22}N_4O_5 \cdot CH_4O_3S$
<b>Purity</b>	greater than or equal to 98% (HPLC)
<b>MW</b>	398.4 . 96.1
<b>CAS Number</b>	59721-29-8
<b>SMILES</b>	<chem>NC(NC1=CC=C(C(OC2=CC=C(CC(OCC(N(C)C)=O)=O)C=C2)=O)C=C1)=N.OS(=O)(C)=O</chem>
<b>Solubility (25°C)</b>	Soluble in water (20mg/ml).
<b>Expiration Date</b>	12 months from date of receipt.

