


## Product Datasheet

# IL2 Recombinant Protein (e.coli Expressed) (orb1227617)

<b>Description</b>	IL2 Recombinant Protein (e.coli Expressed)
<b>Tested Applications</b>	WB
<b>Target</b>	IL2
<b>Preservatives</b>	Lyophilized from 0.22 um filtered solution in PBS, pH 5.0. Normally trehalose is added as protectant before lyophilization
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	Tag free
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for WB. For research use only.
<b>Protein Sequence</b>	Ala 21 - Thr 153
<b>Purity</b>	> 95% as determined by SDS-PAGE.
<b>Source</b>	E. coli
<b>MW</b>	15.4 kDa
<b>Uniprot ID</b>	<b>P60568</b>
<b>NCBI</b>	<b>P60568</b>
<b>Dilution Range</b>	This recombinant protein can be used for WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Human IL-2, Tag Free on SDS-PAGE under r...