


## Product Datasheet

### VEGF R2 Recombinant Protein (orb1227622)

<b>Description</b>	VEGF R2 Recombinant Protein
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	VEGF R2
<b>Preservatives</b>	PBS, pH 7.4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Protein Sequence</b>	Ala 20 - Glu 764
<b>Purity</b>	> 97% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
<b>Source</b>	HEK293 cells
<b>MW</b>	85.2 kDa
<b>Uniprot ID</b>	<a href="#">F7E313</a>
<b>NCBI</b>	<a href="#">F7E313</a>
<b>Dilution Range</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Rhesus macaque VEGF R2, His Tag on SDS-P...