


Product Datasheet

BAFFR Recombinant Protein (orb1227701)

Description	BAFFR Recombinant Protein
Tested Applications	ELISA, WB
Target	BAFFR
Preservatives	Tris with Glycine, Arginine and NaCl, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Ser 10 - Ala 71
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	33.3 kDa
Uniprot ID	Q9D8D0
NCBI	NP_082351.1
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
Human
BAFFR His
Tag at 5
ug/...