


## Product Datasheet

### B7-1 Recombinant Protein (orb1227707)

<b>Description</b>	B7-1 Recombinant Protein
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	B7-1
<b>Preservatives</b>	Tris with Glycine, Arginine and NaCl, pH 7.5
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	Fc Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Protein Sequence</b>	Val 38 - Lys 245
<b>Purity</b>	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
<b>Source</b>	HEK293 cells
<b>MW</b>	50.3 kDa
<b>Uniprot ID</b>	<a href="#">Q00609</a>
<b>NCBI</b>	<a href="#">Q00609</a>
<b>Dilution Range</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Immobilized  
Mouse B7-1,  
Fc Tag at 2  
ug/m...