


Product Datasheet

B7-1 Recombinant Protein (orb1227708)

Description	B7-1 Recombinant Protein
Tested Applications	ELISA, WB
Target	B7-1
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Val 38 - Lys 245
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	25.6 kDa
Uniprot ID	Q00609
NCBI	Q00609
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized Mouse B7-1, His Tag at 2 ug/...