


Product Datasheet

VEGF121 Recombinant Protein (orb1227715)

Description	VEGF121 Recombinant Protein
Tested Applications	ELISA, WB
Target	VEGF121
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Ala 27 - Arg 147
Purity	> 95% as determined by reduced SDS-PAGE. Less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	15 kDa
Uniprot ID	P15692
NCBI	P15692
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
VEGFR1/R2-
Fc at 1
ug/ml can ...