


## Product Datasheet

### PDGF-B Recombinant Protein (orb1227731)

<b>Description</b>	PDGF-B Recombinant Protein
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	PDGF-B
<b>Preservatives</b>	0.085% TFA in 30% ACN
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	His Tag & Avi Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Protein Sequence</b>	Ser 82 - Thr 190
<b>Purity</b>	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
<b>Source</b>	E. coli
<b>MW</b>	15 kDa
<b>Uniprot ID</b>	<a href="#">P01127</a>
<b>NCBI</b>	<a href="#">AAH29822</a>
<b>Dilution Range</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Human PDGF-BB Protein, His Tag, Avi Tag  
...