


Product Datasheet

LILRB3 Recombinant Protein (orb1227743)

Description	LILRB3 Recombinant Protein
Tested Applications	WB
Target	LILRB3
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gly 24 - Glu 443
Purity	> 92% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	48 kDa
Uniprot ID	AAB68668(P288R)
NCBI	AAB68668 , P288
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human LILRB3, His Tag on SDS-PAGE under ...