


Product Datasheet

TGF-beta 1 Recombinant Protein (orb1227747)

Description	TGF-beta 1 Recombinant Protein
Tested Applications	WB
Target	TGF-beta 1
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Leu 30 - Ser 390
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	42.1 kDa (monomer)
Uniprot ID	P01137
NCBI	P01137
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human Latent TGFβ1, His Tag on SDS-PAGE ...