


Product Datasheet

Activin A Recombinant Protein (orb1227748)

Description	Activin A Recombinant Protein
Tested Applications	ELISA, WB
Target	Activin A
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Ser 21 - Ser 426
Purity	> 90% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	13.0 kDa (mature) and 32.0 kDa (pro)
Uniprot ID	P08476
NCBI	AAH07858.1
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
Human
Activin A,
His Tag at ...