


Product Datasheet

LAP Recombinant Protein (orb1227749)

Description	LAP Recombinant Protein
Tested Applications	WB
Target	LAP
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Leu 30 - Arg 278
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	29.8 kDa
Uniprot ID	P01137
NCBI	NP_000651 , C33
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human LAP (TGF-beta 1), His Tag on SDS-P...