


Product Datasheet

IL-37 Recombinant Protein (orb1227750)

Description	IL-37 Recombinant Protein
Tested Applications	WB
Target	IL-37
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Val 46 - Asp 218
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	E. coli
MW	20.3 kDa
Uniprot ID	Q9NZH6
NCBI	AAH20637
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-37b, His Tag on SDS-PAGE under ...