


Product Datasheet

IL-17F Recombinant Protein (orb1227751)

Description	IL-17F Recombinant Protein
Tested Applications	WB
Target	IL-17F
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Arg 31 - Gln 163
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	15.7 kDa
Uniprot ID	Q96PD4
NCBI	AAH70124
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-17F, His Tag on SDS-PAGE under ...