


Product Datasheet

Gremlin Recombinant Protein (orb1227760)

Description	Gremlin Recombinant Protein
Tested Applications	ELISA, WB
Target	Gremlin
Preservatives	50 mM Tris, 150 mM NaCl, pH 7.5 with 1 mM EDTA
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Lys 25 - Asp 184
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	44.2 kDa
Uniprot ID	O60565
NCBI	NP_037504
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized Human BMP-2, Tag Free at 2 u...