


Product Datasheet

BAFFR Recombinant Protein (orb1227796)

Description	BAFFR Recombinant Protein
Tested Applications	ELISA, WB
Target	BAFFR
Preservatives	Tris with Glycine, Arginine and NaCl, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Ser 7 - Ala 71
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	33.2 kDa
Uniprot ID	Q96RJ3
NCBI	Q96RJ3
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
Human
BAFF, His
Tag at 5
ug/...