


Product Datasheet

Angiopoietin-2 Recombinant Protein (orb1227802)

Description	Angiopoietin-2 Recombinant Protein
Tested Applications	WB
Target	Angiopoietin-2
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Lys 275 - Phe 496
Purity	> 95% as determined by reduced SDS-PAGE. Less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	26.3 kDa
Uniprot ID	O15123
NCBI	AAI26201.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human Angiopoietin-2, His Tag on SDS-PAG...