


Product Datasheet

AIMP1 Recombinant Protein (orb1227803)

Description	AIMP1 Recombinant Protein
Tested Applications	WB
Target	AIMP1
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 2 - Lys 312
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	E. coli
MW	35.2 kDa
Uniprot ID	NP_004748.2
NCBI	AAH14051
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human AIMP1, His Tag on SDS-PAGE under r...