


## Product Datasheet

### B7-1 Recombinant Protein (orb1227823)

<b>Description</b>	B7-1 Recombinant Protein
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	B7-1
<b>Preservatives</b>	Tris with Glycine, Arginine and NaCl, pH 7.5
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	Fc Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Protein Sequence</b>	Val 35 - Asn 242
<b>Purity</b>	> 95% as determined by reduced SDS-PAGE. Less than 1.0 EU per ug by the LAL method.
<b>Source</b>	HEK293 cells
<b>MW</b>	50.5 kDa
<b>Uniprot ID</b>	<a href="#">G7NXN7</a>
<b>NCBI</b>	<a href="#">G7NXN7</a>
<b>Dilution Range</b>	This recombinant protein can be used for E, WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Cynomolgus/Rhesus macaque B7-1, Fc Tag o...