


Product Datasheet

Oncostatin M Recombinant Protein (orb1227825)

Description	Oncostatin M Recombinant Protein
Tested Applications	WB
Target	Oncostatin M
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 26 - Arg 252
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	25.8 kDa
Uniprot ID	P13725
NCBI	NP_065391.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human Oncostatin M on SDS-PAGE under red...