


Product Datasheet

G-CSF Recombinant Protein (orb1227826)

Description	G-CSF Recombinant Protein
Tested Applications	ELISA, WB
Target	G-CSF
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant protein can be used for E, WB. For research use only.
Protein Sequence	Thr 31 - Pro 204
Purity	> 95% as determined by SDS-PAGE. Endotoxin level is less than 1.0 EU per ug by the LAL method.
Source	HEK293 cells
MW	18.7 kDa
Uniprot ID	P09919
NCBI	NP_757373.1
Dilution Range	This recombinant protein can be used for E, WB. For research use only.
Expiration Date	6 months from date of receipt.



Human G-CSF, Tag Free on SDS-PAGE under ...