


Product Datasheet

VEGF 121 Recombinant Protein (orb1227837)

Description	VEGF 121 Recombinant Protein
Tested Applications	WB
Target	VEGFA
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 27 - Arg 147
Purity	> 95% as determined by SDS-PAGE of reduced (+) and non-reduced (-) rhVEGF121.
Source	HEK293 cells
MW	14 kDa (monomer)
Uniprot ID	P15692
NCBI	NP_001020541
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human VEGF121, Tag Free on SDS-PAGE unde...