


Product Datasheet

DR4 Recombinant Protein (orb1227850)

Description	DR4 Recombinant Protein
Tested Applications	WB
Target	TNFRSF10A
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 24 - Asn 239
Purity	> 92% as determined by SDS-PAGE.
Source	HEK293 cells
MW	24 kDa
Uniprot ID	O00220
NCBI	NP_003835
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human TRAIL R1, His Tag on SDS-PAGE unde...