


Product Datasheet

Neuropilin 1 Recombinant Protein (orb1227939)

Description	Neuropilin 1 Recombinant Protein
Tested Applications	WB
Target	NRP1
Preservatives	Tris with Glycine, Arginine and NaCl, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Phe 22 - Lys 644
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	96.5 kDa
Uniprot ID	O14786
NCBI	AAH07533
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh NRP1/CD304 Fc Chimera W...