


Product Datasheet

Neuropilin 1 Recombinant Protein (orb1227940)

Description	Neuropilin 1 Recombinant Protein
Tested Applications	WB
Target	NRP1
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Phe 22 - Lys 644
Purity	> 92% as determined by SDS-PAGE.
Source	HEK293 cells
MW	70.7 kDa
Uniprot ID	O14786
NCBI	AAH07533
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human Neuropilin-1, His Tag on SDS-PAGE ...