


Product Datasheet

LEPR Recombinant Protein (orb1227967)

Description	LEPR Recombinant Protein
Tested Applications	WB
Target	LEPR
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Phe 22 - Asp 839
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	121 kDa
Uniprot ID	P48357
NCBI	NP_002294.2
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh LEPR Fc Chimera was det...