


## Product Datasheet

### CD360 Recombinant Protein (orb1227976)

<b>Description</b>	CD360 Recombinant Protein
<b>Tested Applications</b>	WB
<b>Target</b>	IL21R
<b>Preservatives</b>	50 mM tris, 100 mM glycine, pH 7.5
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	C-Fc Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for WB. For research use only.
<b>Protein Sequence</b>	Cys 20 - Pro 236
<b>Purity</b>	> 92% as determined by SDS-PAGE.
<b>Source</b>	HEK293 cells
<b>MW</b>	51.6 kDa
<b>Uniprot ID</b>	<a href="#">Q9HBE5</a>
<b>NCBI</b>	<a href="#">AAH04348</a>
<b>Dilution Range</b>	This recombinant protein can be used for WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



The purity of rh IL21R/CD360 Fc Chimera ...