


Product Datasheet

IL1R3 Recombinant Protein (orb1227982)

Description	IL1R3 Recombinant Protein
Tested Applications	WB
Target	IL1RAP
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ser 21 - Glu 359
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	40 kDa
Uniprot ID	Q9NPH3
NCBI	Q9NPH3
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh IL1RAP/IL-1RACP/IL1R3 W...