


Product Datasheet

Interleukin 29 Recombinant Protein (HEK293 Expressed) (orb1227992)

Description	Interleukin 29 Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	IFNL1
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Met 1 - Thr 200
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	21.1 kDa
Uniprot ID	Q8IU54
NCBI	NP_742152.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh IL29 was determined by ...