


Product Datasheet

Interleukin 6 Recombinant Protein (HEK293 Expressed) (orb1227994)

Description	Interleukin 6 Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	IL6
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Val 30 - Met 212
Purity	> 98% as determined by SDS-PAGE.
Source	HEK293 cells
MW	20.8 kDa
Uniprot ID	P05231
NCBI	NP_000591
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-6, Tag Free on SDS-PAGE under r...