


Product Datasheet

Interleukin 4 Recombinant Protein (HEK293 Expressed) (orb1227996)

Description	Interleukin 4 Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	IL4
Preservatives	20 mM Na ₂ HPO ₄ , 200 mM NaCl, pH 7.0
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	His 25 - Ser 153
Purity	> 90% as determined by SDS-PAGE.
Source	HEK293 cells
MW	14.9 kDa
Uniprot ID	P05112
NCBI	AAH67514
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-4, Tag Free on SDS-PAGE under r...