


Product Datasheet

IL 1 R2 Recombinant Protein (orb1227998)

Description	IL 1 R2 Recombinant Protein
Tested Applications	WB
Target	IL1R2
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Thr 14 - Glu 343
Purity	> 98% as determined by SDS-PAGE.
Source	HEK293 cells
MW	38.6 kDa
Uniprot ID	P27930
NCBI	NP_004624.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-1 RII, His Tag on SDS-PAGE unde...