


Product Datasheet

Interferon gamma Recombinant Protein (HEK293 Expressed) (orb1228013)

Description	Interferon gamma Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	IFNG
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gln 24 - Gln 166
Purity	> 90% as determined by SDS-PAGE.
Source	HEK293 cells
MW	16.8 kDa
Uniprot ID	P01579
NCBI	AAH70256
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IFN-gamma, Tag Free on SDS-PAGE un...