


Product Datasheet

IFNAR 1 Recombinant Protein (orb1228017)

Description	IFNAR 1 Recombinant Protein
Tested Applications	WB
Target	IFNAR1
Preservatives	Tris with Glycine, Arginine and NaCl, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Lys 28 - Lys 436
Purity	85% (rh IFNAR1 /IFNAR Fc Chimera) and 13% (Human IgG1 Fc fragment) as determined by SDS-PAGE.
Source	HEK293 cells
MW	73.8 kDa
Uniprot ID	P17181
NCBI	NP_000620
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh IFNAR1/IFNAR Fc Chimera...