


Product Datasheet

Influenza A Hemagglutinin Recombinant Protein (orb1228027)

Description	Influenza A Hemagglutinin Recombinant Protein
Tested Applications	WB
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Asp 19 - Gly 338, Leu 341 - Asp 523
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	58.6 kDa
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of HA (H7N9)-Shanghai was det...