


Product Datasheet

TNFRSF 18 Recombinant Protein (orb1228039)

Description	TNFRSF 18 Recombinant Protein
Tested Applications	WB
Target	TNFSF18
Preservatives	50 mM tris, 100 mM glycine, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	N-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gln 50 - Ser 177
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	46.2 kDa
Uniprot ID	Q9UNG2
NCBI	AAH69319.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human
GITR
Ligand, Fc
Tag on SDS-
PAGE un...