


Product Datasheet

Flt3 Ligand Recombinant Protein (HEK293 Expressed) (orb1228053)

Description	Flt3 Ligand Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	FLT3LG
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Thr 27 - Pro 185
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	18.9 kDa
Uniprot ID	P49771
NCBI	AAI26294
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human Flt-3 Ligand, Tag Free on SDS-PAGE...