


Product Datasheet

E Cadherin Recombinant Protein (orb1228079)

Description	E Cadherin Recombinant Protein
Tested Applications	WB
Target	CDH1
Preservatives	50 mM tris, 100 mM glycine, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gln 23 - Pro 621
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	92.4 kDa
Uniprot ID	P12830
NCBI	AAI41839
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rhECAD Fc Chimera was dete...