


Product Datasheet

CXCR4 Recombinant Protein (orb1228093)

Description	CXCR4 Recombinant Protein
Tested Applications	WB
Target	CXCR4
Preservatives	50 mM tris, 100 mM glycine, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	N-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Met 1 - Ser 46
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	32.8 kDa
Uniprot ID	P61073
NCBI	AAH20968.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh CXCR4 Fc Chimera was de...