


## Product Datasheet

### IL 2 R beta Recombinant Protein (orb1228156)

<b>Description</b>	IL 2 R beta Recombinant Protein
<b>Tested Applications</b>	WB
<b>Target</b>	IL2RB
<b>Preservatives</b>	PBS, pH 7.4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
<b>Tag</b>	His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for WB. For research use only.
<b>Protein Sequence</b>	Ala 27 - Asp 239
<b>Purity</b>	> 90% as determined by SDS-PAGE.
<b>Source</b>	HEK293 cells
<b>MW</b>	25.4 kDa
<b>Uniprot ID</b>	<a href="#">P14784</a>
<b>NCBI</b>	<a href="#">NP_000869.1</a>
<b>Dilution Range</b>	This recombinant protein can be used for WB. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.



Human IL-2 R beta, His Tag (SPR verified...)