


Product Datasheet

CD80 Recombinant Protein (orb1228163)

Description	CD80 Recombinant Protein
Tested Applications	WB
Target	CD80
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 37 - Ser 160
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	25.8 kDa
Uniprot ID	G7MKE2
NCBI	G7MKE2
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Cynomolgus/Rhesus macaque B7-1, His Tag ...