


Product Datasheet

IL-33 Recombinant Protein (orb1228217)

Description	IL-33 Recombinant Protein
Tested Applications	ELISA, WB
Target	IL33
Preservatives	100mM sodium phosphate, 500mM NaCl, 25 mM imidazole, 10% glycerol
Form/Appearance	Liquid
Storage	Store in working aliquots at -70°C. Avoid freeze/thaw cycles. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	Fusion Partner: C-terminal His-tag
Note	For research use only
Application notes	This recombinant protein can be used for WB and ELISA. For research use only.
Protein Sequence	aa 112 - 270
Purity	~95%
Source	E. coli
MW	22 kDa (Calculated)
Uniprot ID	O95760
NCBI	NP_254274
Dilution Range	This recombinant protein can be used for WB and ELISA. For research use only.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of recombinant IL-33 f...