

## Product Datasheet

# SARS-CoV-2 (COVID-19) Helicase Recombinant Protein (orb1228226)

<b>Description</b>	SARS-CoV-2 (COVID-19) Helicase Recombinant Protein
<b>Target</b>	ORF1ab
<b>Preservatives</b>	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.5.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store below -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	6His tag at the C-terminus
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Ala1-Gln601
<b>Purity</b>	Greater than 80% as determined by reducing SDS-PAGE.
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<a href="#">P0DTC1</a>
<b>NCBI</b>	<a href="#">YP_009725297.1</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.