


Product Datasheet

Rhesus Macaque ACE-2 Recombinant Protein (orb1228245)

Description	Rhesus Macaque ACE-2 Recombinant Protein
Target	Ace2
Preservatives	Supplied as a 0.2 um filtered solution of 25mM Tris-HCl, 150mM NaCl, 1mM ZnCl ₂ , pH 7.5.
Form/Appearance	Liquid
Storage	Store below -20°C, stable for 6 months after receipt.
Tag	C-10 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Gln18-Val739
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Mammalian Cells
MW	85.1kDa
Uniprot ID	ACI04564.1
NCBI	ACI04564.1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
Rhesus
Macaque
ACE-2-His
(Ca...)